

Explore Genetic Technology to Improve the Quality of Life

Pluslife CRQ Card

One Nucleic Acid Test Card Four Canine Pathogen Identification



CPIV

CAV-2



1 1 1

pluslife




Features

Pluslife Mini Dock is a rapid molecular diagnostic system that can be used to detect canine distemper virus (CDV), canine parainfluenza virus (CPIV), canine adenovirus type 2 (CAV-2), *Brodetella bronchiseptica* (Bb) in one single test. Our self-developed isothermal amplification technology, RHAM, provides excellent results comparable to qPCR tests. The product is compact and portable, with easy-to-operate procedure to obtain results within 35 minutes.


FAST TAT within 35 mins Image: Construction of the structure of the structure

Clinical Proven



| РРА |
|---------------------|
| CDV: 100% |
| CAV-2: 94.7% |
| CPIV: 100% |
| Bb: 100% |

Test Procedure



Result Output



Application



Pet Hospital/Clinic




Pet Shop



Center for Animal Disease Control

Specification





| Pluslife Mini Dock | | |
|------------------------|---------------------------------|--|
| Dimension | 101 x 91 x 65 mm | |
| Weight | 210 g | |
| Connection | Mobile phone or computer app | |
| Operate temperature | 15°C–30°C | |
| Operate humidity | ≤80% | |

| Pluslife CRQ Card | | |
|----------------------------------|--|--|
| CDV, CPIV, CAV-2, Bb | | |
| Eye/nasal/ oropharyngeal swab | | |
| 1000 copies/mL | | |
| 35 mins | | |
| 13 months at room temperature | | |
| | | |

Guangzhou Pluslife Biotech Co.,Ltd.



PLCRQ-INT-202303020